

ProCare Medical

healthcare solutions²



Lopital *Elexo* Shower Chair

More and more care facilities are turning to showering for resident's personal hygiene. Showering a resident in a standard chair will cause the caregiver to work in unnatural postures and suffer static back strain. The best way to overcome this situation is to introduce the right equipment that significantly increases the percentage of time that the caregiver spends working in a healthy position. Lopital offers relief with their extensive line of shower chairs and trolleys.

The ***Elexo*** shower chair provides the resident maximum safety and comfort and a single caregiver can take care of dressing and undressing, showering, toileting, foot care, and hair washing in an ergonomically sound working condition.



Over toilet seat



Easy-to-use hand control



17 ° tilt adjustment



Features

- ◆ Electrically operated height-adjustable settings.
- ◆ Easy-to-use hand control with suspension hook.
- ◆ Tilt adjustment to a maximum of 17 °.
- ◆ Fold-out footrests with calf protection.
- ◆ Fold-away and removable leg supports.
- ◆ Ergonomically designed backrest and seat.
- ◆ Foldable padded armrests with front lock for lateral transfer.
- ◆ Suitable for use in the shower and can be moved into position over toilet.
- ◆ Guide for plastic commode pan.
- ◆ Easy to recharge with battery charger provided.
- ◆ Wide wheelbase for stability.
- ◆ Large diameter swivel castors with brakes.
- ◆ Corrosion resistant stainless steel, galvanized, and polymer parts.
- ◆ Extremely light aluminum frame.

Technical specifications

- ◆ Weight capacity: 400 lbs.
- ◆ Weight: 110 lbs.
- ◆ Tilt adjustment: 0 - 17 ° maximum
- ◆ Overall dimensions (frame): 32.7" L x 26.8" W
- ◆ Adjustable sitting height: 19.6" minimum / 41.3" maximum
- ◆ Seat dimensions: 18.9" W / 18.1" Depth
- ◆ Castor diameter: 5"
- ◆ Electric controls: 24V
- ◆ Battery: 12V 7Ah,
- ◆ Battery charger: 220V 50/60Hz
- ◆ Number of lifts (fully charged battery): 40 up-down movements

