

ProCare Medical

healthcare solutions²

Handling belts are used as an alternative to holding the users clothing or body while standing or walking. Handling belts can be used when assisting with standing from a chair or bed and for general walking transfers. The user must be able to maintain and support their own weight when using this product. It must not be used as a standing or lifting aid.



The **Deluxe Patient Belt** is manufactured in a quilted polyester fabric with a non-slip inner lining and a click fit buckle with an extra wide back to improve the wearer's comfort. The Deluxe Patient belt features a variety of grab handles for ease of use.



The **Mini Deluxe Patient Belt** is constructed of the same quilted polyester fabric with the same non-slip inner lining and click fit buckle as the Deluxe Patient Belt, but is of a uniform comfortable 4 inch width for the entire length of the belt

Sylvia Patient Belt Sizing

Standard for waists 26" to 36" (Not available in Deluxe)*

Medium-Large for waist size: 39"-47" *

XL now available! for waist size: 45" - 51" *

- * For maximum comfort, to allow the belt to overlap on itself slightly, use 26" maximum waist measurement for the Standard, use 38" maximum waist measurement for the Medium/Large and 44" maximum waist measurement for the Large/ XL