



RR7-II Atlantic



RG9 Victoria



RM3-I Superior

ProCare Medical

healthcare solutions ²



RL9 Bath Lift



RT4 Chesapeake



RS8 Geneva

Bathing Systems for Healthcare Facilities

Rane Bathing Systems has been manufacturing quality bathing systems for care facilities since 1999. Bathing promotes a feeling of well-being and is mentally and physically beneficial, which is why Rane is dedicated to producing innovative bathing systems. Rane bathing systems are designed to be functional and affordable, yet will give your bathing area an upscale look. Featuring numerous standard features and options, Rane bathing systems will enhance the bathing experience for the resident.

The Rane Bathing Systems product line provides a unique combination of style, ease of use, comfort, and low cost of ownership.

- ◆ **100% Made in the USA.** Manufactured in Sparta, TN from the highest quality components and materials.
- ◆ **A “Best in Class” Warranty** resulting in low cost of ownership.
- ◆ **Full set of Standard Features**, including the CleanRane™ Air Spa System, thermostatic mixing valve, handheld shower wand, built-in disinfection system, easy-to-use control panel, water temperature display, and 3/4” stainless steel mesh supply lines.
- ◆ **Nationwide Distribution and Customer Service.**



All Rane bathing systems feature this easy-to-use control panel which includes the CleanRane™ Air Spa control, thermostatic mixing valve, tub-fill valve, shower valve, water temperature display, and handheld shower wand.



The air jets on the CleanRane™ Air Spa System feature spring loaded check valves for back flow prevention, eliminating the extra step to disinfect the air lines.



100% MADE IN THE USA



RR7-II Atlantic Reclining Bathing System

- ◆ At only 24 gallons to fill the foot well, water usage and the associated costs are greatly reduced.
- ◆ Unique tilting function, which allows patients with challenging conditions to recline in a contoured seat and become immersed in warm water.
- ◆ 3 plumbing configurations to choose from-side, front, or island configuration.
- ◆ Optional Auto-fill is a time-saving feature that enables the caregiver to simply press a button and pre-fill the tub to a pre-set level.



RG9 Victoria Side-Entry Bathing System

- ◆ Designed specifically for use with the RL9 Victoria Bath Lift, which provides seamless transfers from bed-to-bath in a comfortable seated position.
- ◆ The low threshold and wide side-entry door allows easy access for residents with various degrees of mobility.
- ◆ The wide built-in contoured seat allows the larger resident to sit comfortably and safely in a slightly reclined position.
- ◆ 600 lb. weight capacity.



RS8 Geneva Supine Bathing System

- ◆ Designed to bathe highly dependent residents with challenging conditions in a comfortable supine position.
- ◆ Available in either Height Adjustable or Fixed Height.
- ◆ Compatible with a wide range of transport and lifting devices.
- ◆ Optional Auto-fill is a time-saving feature that enables the caregiver to simply press a button and pre-fill the tub to a pre-set level.



RT4 Chesapeake Built-In Bathing System

- ◆ Designed to be recessed into a standard bathtub opening.
- ◆ Ideal for care facilities where residents can bathe unassisted.
- ◆ Wall-mounted modular control panel that can be positioned for convenient operation by either the resident or caregiver.
- ◆ Optional wall surround made from the same high-quality fiberglass and gel coat as the RT4 Chesapeake.



RM3-I Superior Built-In or Freestanding Bathing System

- ◆ Designed to be a freestanding, or built-in, side access bathing system.
- ◆ Measuring just 30" wide by 50" in length, the RM3-I Superior is designed for care facilities with small bathing areas or limited space.
- ◆ The unique door is designed to clear commodes with ease.
- ◆ Optional Wall-Mounted Disinfection Cabinet available.

Rane's 'Best in Class' Warranty

- ◆ 5 years for external finishes and composite shell
- ◆ 3 years on all moving and stationary parts
- ◆ Lifetime warranty on the door seal

UL/CSA Certifications: ANSI Z21.20 / UL 1795 | CSA C22.2 No. 218.2 | ANSI Z124.1.2 (2005): CAN / CSA B45

ProCare Medical | 1650 Manheim Pike, Suite 102, Lancaster, PA 17601 | 855-528-0421 | Fax: 877-528-0421

www.procare-medical.com