



**BRODA**<sup>®</sup>  
*Comfort is Essential*

**SYNTHESIS**  
TILT RECLINER V4 ACC

## ENHANCED INFECTION CONTROL FOR MULTI-PATIENT CARE

Designed for acute care and hospital settings, the V4 ACC combines enhanced infection control features, with industry leading functionality, making it the unrivaled solution for your next clinical application.



- 1 INFINITELY ADJUSTABLE SEAT TILT
- 2 INFINITELY ADJUSTABLE BACK RECLINE
- 3 HEIGHT ADJUSTABLE ARMRESTS
- 4 SWING AWAY REMOVABLE ARMS
- 5 ADJUSTABLE WINGS FOR UPPER LATERAL SUPPORT
- 6 INDEPENDENT ELEVATING, ARTICULATING, LENGTH ADJUSTABLE LEGREST
- 7 ABS FLIP DOWN FOOTREST
- 8 MULTI-DIRECTIONAL LOCKING REAR CASTERS

### PATIENT BENEFITS



Provides enhanced infection control with fully wrapped seat, back and legrest covers.



Allows patients to be safely transported throughout the facility.



Facilitates ease of transfers, accommodating safe patient handling and no-lift programs.



Wings provide height and width adjustability for individualized upper lateral support.

### CAREGIVER AND FACILITY BENEFITS

- Design provides safety and stability helping to reduce risk of injury
- Fully wrapped covers provide ease of cleaning in multi-patient care environments
- Posterior tilt enhances pressure redistribution aiding in the maintenance of skin integrity
- Tilt and recline provide ease of use for frequent re-positioning, helping improve patient outcomes
- Enhanced positioning and innovative features allow for fewer transfers, promoting safe patient handling
- Quality construction provides superior durability, reducing equipment maintenance and replacement costs

## COMFORT TENSION SEATING® PROVIDES

- Long term seating comfort
- Conforming support for proper pelvic positioning
- Improved air flow, decreasing heat and moisture build up
- Pressure redistribution that aids in maintaining skin integrity
- Accommodation for conditions such as kyphosis and lordosis
- Decreased pressure, improved blood flow and oxygenation of the tissue



CTS Strapping flexes for individualized fit & superior comfort.


## SPECIFICATIONS

Frame	16 Gauge tubular steel frame
Seat & Backrest	BRODA Comfort Tension Seating®
Degrees of Tilt	40°
Degrees of Recline	90° (Trendelenburg capable)
Standard Seat Width	16", 18", 20" (40cm, 45cm, 50cm)
Custom Seat Width	22", 24", 26" (55cm, 60cm, 65cm)
Standard Seat Depth	17", 18.5" (43cm, 46cm)
Custom Seat Depth	20" (50cm)
Standard Seat Height	20.5" (51cm)
Custom Seat Height	19", 22", 23.5" (48cm, 55cm, 59cm)
Front Casters	5" Front total lock (12.5)
Rear Caster	6" Multi-directional locking
Mag Wheel Upgrade	20" or 22" Rear mount (50cm or 55cm), 14" rear or mid mount (35cm)
Weight Capacity	350lbs. (159kg)
Standard Padding	Headrest cushion, shoulder bolsters, lower side pads, fully wrapped seat and back pad, calf pad, sole pads

## COMMON ACCESSORIES


- Mag wheel package
- Aluminum front casters
- Solid or clear tray
- Half tray
- Padded tray cover
- Padded armrests
- O2 Holder
- IV Pole

## UPHOLSTERY OPTIONS



\* Out of respect for their service camo options available for veterans only

## FRAME OPTIONS





HEALTHCARE DIVISION

[info@brodaseating.com](mailto:info@brodaseating.com)

p. 1-844-552-7632

[Brodaseating.com](http://Brodaseating.com)

©Broda USA, INC.

8/16