

# Medcare Car Extractor



The Medcare Car Extractor is designed to safely lift and remove a patient from nearly any vehicle ranging from a compact automobile to fire or emergency van. It is ideal for hospital emergency departments as this heavy-duty lift enables speedy transfers of patients onto gurneys, stretchers, or wheelchairs. The Car Extractor can also insert patients into vehicles effectively.

Built light and strong with a 30" reach and a 40" lifting range, the Car Extractor can gain deep access into tight vehicle interior spaces and lift patients with limited mobility or little weight bearing ability—up to 475 lbs. capacity.

Item Number:

450000: Car Extractor 475 -  
lifts up to 475 lbs.

450045: Low Pro Car Extractor 475 -  
lifts up to 475 lbs.



“ Perfect for reducing injuries related to vehicle transfers, the Car Extractor is ideal for transferring patients in and out of most vehicles. ”



# Car Extractor (cont.)



### Hanger Bar Stopper Option

The Medicare Hanger Bar Stopper prevents sling loops from sliding off of the hanger bar hooks.

Item Number: 400611



Hand Control



Battery Charger

Battery

## Technical Specifications

Overall Length:	47"	Boom Arm Reach:	30"
Outside Width (legs closed):	27"	Lifting Range:	40"
Inside Width (legs open):	36"	Power Source:	Rechargeable 24 Volt, 5 Amp Battery
Height of Wheel Base:	4.5"		

### Product Highlights:

- Dual battery design allows stand to be in constant use
- 360° rotation for easy patient positioning
- Straight boom for deep access into vehicle
- Hand control allows caregiver to be positioned next to patient during transfer
- Compact, padded 2-point hanger bar
- To be used with two Medicare Stand Belts and two personnel



### ProCare Medical

1650 Manheim Pike,  
Ste 102  
Lancaster, PA 17601

Ph: 855-528-0421  
Fx: 877-528-0421

[www.procare-medical.com](http://www.procare-medical.com)