

Comparison of Volker to Medmizer Bed



Item	Volker	Med-Mizer
Manufacturing	Germany, now part of Hillrom	Batesville, Indiana
Controls	Handcontrol	Handcontrol, or integrated to assist bars
Section design	4 section, butterfly deck	4 section, with clinical contour
Low heights	13"	9"
Widths	34" or 38"	35" standard 35-39-42" expanding feature 35-42-48" bariatric unit, expanding feature
Mattresses	European sizes	Fits standard size
Lead time	Was 8-12 weeks	1 week
Wall Guard	Optional to protect wall from headboard	Included on head board as wheel bumpers
Pushing bed	Could only be maneuvered at lowest position	Can push at any height
Storage/Service	Flat	Vertical
ReTractability	NA	Included
Side Rails	Included, telescoping and fold up from bed frame	Assist bar or ¼" length rail

Ask for the complete Med-Mizer Catalog and Brochure for complete details.

Click for video of clinical feature of bed [here](#)