

# ProCare Medical

*healthcare solutions<sup>2</sup>*

## Showering in a supine position

For some people with a disability, it is more pleasant or better to be washed and taken care of in a lying position. The Sirocco allows this. This enables an immobile person to be cared for in comfort.

The correct shower table will ensure that the caregiver has greater freedom of movement, can work closer to the care recipient and the work causes less strain. It also enables the care recipient to play a more active role in his or her personal care.

Depending on the space that is available and the needs, you can choose from several wallmounted stretchers. There is also a choice between electric and fixed height shower tables.

For more information  
visit our website

[www.procare-medical.com](http://www.procare-medical.com)



6100 2400

# Sirocco Shower stretcher

## Technical specifications

| Product no. | Product Description                          | Outside Length | Outside Width | Inside length | Inside width | Lying surface height adjustment | Max. load kg | Max. load lbs | Product weight in kg | Product weight in lbs |
|-------------|--|----------------|---------------|---------------|--------------|---------------------------------|--------------|---------------|----------------------|-----------------------|
| 6100 2400   | Sirocco Shower stretcher wall-mounted, left  | 78.75          | 32            | 74.80         | 23.6         | 17 - 40.5                       | 150          | 330           | 80                   | 180                   |
| 6100 2401   | Sirocco Shower stretcher wall-mounted, right | 78.75          | 31            | 74.80         | 23.6         | 17 - 40.5                       | 150          | 330           | 80                   | 180                   |

All measurements in this matrix are in inches.



- Wall-mounted.
- Foldable.
- Shock-free electric adjustment from 40 to 100 cm.
- Only suitable for bathrooms with a minimum height of 230 cm.



Folding down side rail.



Removable battery.

# Amfora Shower stretcher



6300 4030 / 6300 4035  
**Amfora Shower stretcher wall mounted, hinged without collection tray**



6300 4040 / 6300 4045  
**Amfora Shower stretcher wall mounted, hinged with collection tray**

## Technical specifications

| Product no. | Product Description                                      | Outside Length | Outside Width | Inside length | Inside width | Lying surface height adjustment | Max. load kg | Max. load lbs | Product weight in kg | Product weight in lbs |
|-------------|--|----------------|---------------|---------------|--------------|---------------------------------|--------------|---------------|----------------------|-----------------------|
| 6300 4030   | Amfora Shower stretcher wall 160 without collection tray | 63             | 33.5          | 59            | 23           | 33.5                            | 135          | 300           | 18                   | 40                    |
| 6300 4035   | Amfora Shower stretcher wall 190 without collection tray | 75             | 33.5          | 71            | 23           | 33.5                            | 135          | 300           | 18                   | 40                    |
| 6300 4040   | Amfora Shower stretcher wall 160 with collection tray    | 63             | 33.5          | 59            | 23           | 33.5                            | 135          | 300           | 20                   | 44                    |
| 6300 4045   | Amfora Shower stretcher wall 190 with collection tray    | 75             | 33.5          | 71            | 23           | 33.5                            | 135          | 300           | 20                   | 44                    |

All measurements in this matrix are in inches.