

Vienna Dressing Tables



6500 4021 / 6500 4026
Vienna Dressing Table
 electric height adjustable




6500 4020 / 6500 4025
Vienna Dressing Table
 hydraulic height adjustable


Technical specifications 6500 4020

- Max. load 150 kg / 330 lbs
- Product weight 56 kg / 124.0 lbs
- Length 160 cm / 63.0 inch
- Width 82 cm / 32.3 inch
- Height (adjustable) 56 to 106 cm / 22.0 to 41.7 inch

Technical specifications 6500 4025

- Max. load 150 kg / 330 lbs
- Product weight 59 kg / 130.1 lbs
- Length 190 cm / 74.8 inch
- Width 82 cm / 32.3 inch
- Height (adjustable) 56 to 106 cm / 22 to 41.7 inch





6500 5015 / 6500 5020

**Vienna Dressing Table
electric height adjustment**

- Wall mounted.
- Synthetic blue leather upholstery.
- Foldable, gas spring assisted.
- Suitable for bathrooms with a minimum height of 230 cm / 90.6 inch
- Synthetic screen guard.
- Adjustable backrest.
- Side rail.
- **Optional:** extra battery pack.

Technical specifications 6500 5015

• Max. load	150 kg / 330 lbs	• Max. load	150 kg / 330 lbs
• Product weight	75 kg / 165 lbs	• Product weight	78 kg / 172 lbs
• Length	160 cm / 63.0 inch	• Length	190 cm / 74.8 inch
• Width	91 cm / 35.8 inch	• Width	91 cm / 35.8 inch
• Height (adjustable)	49 to 109 cm / 19.3 to 42.9 inch	• Height (adjustable)	49 to 109 cm / 19.3 to 42.9 inch
• High / low drive	Actuator 24V DC	• High / low drive	Actuator 24V DC
• Controlbox	24V DC, max. 250W	• Controlbox	24V DC, max. 250W
• Battery	24V, max. 250W	• Battery	24V, max. 250W
• Charger	230V AC, connection value 100W, wall mounted	• Charger	230V AC, connection value 100W, wall mounted
• Waterproofness	IPX4	• Waterproofness	IPX4