

RAIZER[®] M

Introducing the new, all manual, lifting chair!

Based on market demands, Liftup has developed a new Raizer, specifically for the private consumer market.

Independent living without the need of caretaking from outside

Always Ready

The Raizer M is not depending on electrical power and therefore always ready to be used.

Ease in use

The Raizer M is intuitive and easy to assemble and operate. Raizer M does not require any introduction or training other than the documentation supplied with the product and access to online information.

Weight

The Raizer M has a total weight of 27.3 lb. Heaviest item is the seat, weighing 15.9 lb. This makes the product easy to carry to the area where a person has fallen.

Comfortable lifting

The lifting process is very comfortable compared to other solutions available on the market. This has been proven with Raizer and Raizer II.

Hygiene

All surface materials and structures have been developed to be easy to clean and disinfect.

ProCare Medical
healthcare solutions²

759 Flory Mill Road, Lancaster PA 17601
Ph: 855-528-0421
Fx: 877-528-0421
Web: www.procare-medical.com



What's included

- Lifting Chair
- Crank
- Safety belt
- Headrest

Accessories

- Carry bags
- Trolley
- Stair wheels
- Hygiene cover
- Safety belt
- Wall bracket

Operation

Hand operated by one person

Segment

Targeted towards the private consumer market

Technical specifications

- Lifting capacity: 330 lb.
- Duty cycle: 3-5 min.

