

”More independence
with your personal care”



ROPOX TOILET SUPPORT ARMS

GET SOMETHING TO LEAN ON

Toilet support arms are a great aid if you need to get on and off the toilet safely. Their ergonomic shapes provide a comfortable grip, and the cleaning-friendly surfaces optimize hygiene.

Because they are designed with two heights and an anti-slip surface, you have a firm grip from both standing and seated position. If you need a little extra support for getting up and sitting down, you can mount the standing/sitting attachments.

The support arms can be adjusted in height and at different angles to suit the individual user's needs.

If you need extra room in the bathroom, the support arms can be folded up to vertical position.






ROPOX TOILET SUPPORT ARMS

Height adjustment: 4.7"

Max. user weight: Tested acc. to DS/ISO 17966:2016
(see max. user weight below)




VARIANT/MODEL


WAVE - 35.5"		
		
418 lbs Item no.: 40-40115	860 lbs Item no.: 40-40125	848 lbs Item no.: 40-40145

WAVE - 30"

562 lbs Item no.: 40-40151


Can also be used for Ropox Toilet Lifter

STRAIGHT - 35.5"		
		
418 lbs Item no.: 40-40110	860 lbs Item no.: 40-40120	848 lbs Item no.: 40-40140


STRAIGHT - 30"

Item no.: 40-40150/40-40152

ACCESSORIES

Height adjustment
(adjusted with tool)
Item no.: 40-40906



Height adjustment
(adjusted without tool)
Item no.: 40-40907



Toilet roll holder
Item no.: 40-40910



Standing/Sitting attachment
Item no.: 40-40925

