



Creating balance in healthcare



Tango XXL Shower-toilet chair



5100 5800

Tango XXL Shower-toilet chair



5100 5800

Tango XXL Shower-toilet chair

The Lopital Tango XXL is a shower-toilet chair providing obese clients care in a responsible manner. The Tango provides optimum comfort for the caregiver as well as the client.

The Tango XXL can be used as a shower chair as well as a toilet chair that can be wheeled over the toilet. The extra seat width, depth and the attractive design makes the Tango XXL chair a nice looking and user-friendly product caring for obese clients.

The ergonomic shaped backrest and seat provide comfort and security. The extended armrests can fold away completely. The large comfortable footrests including the calf support can be swiveled or removed allowing the client to easily sit down or stand up.

- Maximum load 360 kg / 800 lbs.
- Mobile
- 3 plastic castors fitted with brakes
- 1 direction wheel
- Rail for bedpan or toilet bucket
- Extended and fully foldable armrests
- Two-piece seat opened at two sides
- Swing-out and removable height adjustable footrests

Options

- Product no. 6950 1104 Plastic bedpan and lid
- Product no. 6975 8044 Safety belt
- Product no. 6975 8045 Four point safety belt






Suitable for people with a mild to moderate disability



Two parted seat with front and back openings. Easy removable snaplock seat.



Independent, foldable armrests.



Direction wheel.




Swing-out and removable height adjustable footrests.



Technical specifications

- | | | | |
|------------------|-------------------|---------------------------|--------------------|
| • Max. load | 360 kg / 800 lbs | • Seat height | 51 cm / 20.1 inch |
| • Product weight | 31 kg / 68.3 lbs | • Total length base frame | 97 cm / 31.1 inch |
| • Seat width | 70 cm / 27.6 inch | • Total width base frame | 83 cm / 32.7 inch |
| • Seat depth | 54 cm / 21.3 | • Turning radius | 115 cm / 45.3 inch |



759 Flory Mill Road, Lancaster PA 17601
 Toll Free: 855-528-0421 | www.ProCare-innovations.com

