


RETRACTABED® Series

COMFORT WIDE™ (Expandable)

PUSH-BUTTON SLIDE TECHNOLOGY
WITH LOCKING POSITIONS 35-39-42
INCHES WIDE



RETRACTABED®
Model

EVAC POSITION™

PASS THROUGH DOORS NO MATTER
THE ASSIST DEVICE ATTACHED

- **Clinical-Contour®** – exclusive feature that reduces skin shear and the repositioning of residents
- **Evac Position™** – extra position to make the bed thinner so if a side bar is attached it can still pass through doors less than 42 inches wide (**expandable only**)
- **Comfort Wide™** – bed expands to 42 inches wide without the use of tools for more comfort and flexibility (**expandable only**)
- **Staff Lockout Control** – standard feature built into pendant to prevent unnecessary operation
- **Wall Bumpers** – soft non-marking bumpers come standard to protect wall from damage
- **Retracting Bed Frame** – reduces potential falls and increases staff productivity
- **Embedded Bed Controls** – optional fixed controls in assist bar or head rail reduces replacement handsets and increases staff productivity
- **SureLock** – optional one step lock for improved ease of use while locking casters



SAFETY IN CAREGIVING FOR HEALTHCARE

759 Flory Mill Road, Lancaster PA 17601
Ph: 855-528-0421
Fx: 877-528-0421
www.procare-medical.com

MED MIZER®

RetractaBed® Model SPECIFICATIONS

Deck Height Low	9"
Deck Height High	25"
Expandable Sleep Surface.....	35"w, 39"w, 42"w x 80"l
Fixed Sleep Surface	35"w x 80"l
Maximum Safe Working Capacity	500 lbs
Product Weight	225 lbs

LIMITED WARRANTY

15 years on bed frame mechanics	1 year on wood, casters & handsets
5 years on drive system	Lifetime warranty on all welds



ENDBOARD OPTIONS



ACCESSORIES



PROUDLY MADE IN THE U.S.A.

(855) 528-0421
ProCare Medical

Approvals: IEC60601-1:2005 2nd edition, ANSI/AAMI ES60601-1:2005 2nd edition, and CAN/CSA-22.2 No 60601-1:2008



is at the forefront of innovating medical technology for post acute care. We are focused on partnering with you to help address healthcare challenges in the following key areas—

