

Powered Surfaces (I-MersAir Mattress)

Med-Mizer Powered Surface

- Low air loss and alternating pressure with 5 selectable cycle time and CLP
- 8 settings of firmness/softness for comfort
- Auto Firm mode assist with transfers and repositioning
- 20 8" clear "soft" air bladders for added comfort and infection control
- Cell-on-cell technology
- Constant low pressure mode (CLP) for advanced immersion and envelopment
- Auto detection sets appropriate pressure output for patients weight
- Two-way stretch cover with welded seams (anti-fungi, antibacterial, vapor permeable, water resistant)
- 36" (400 lb capacity) and 42" (550 lb capacity)



MED MIZER®