

Midline

Positioning Wheelchair

POSITIONING
BRODA[®]



Full tilt and recline functionality provides superior positioning




Adjustable seat surface height adapts the Midline to support users at multiple care levels



Tilt-in-space enhances comfort

Patient, Facility, & Caregiver Benefits



- Increases quality of life and engagement with surroundings
- Allows for easy maneuverability and safe patient adjustments
- Reduces risk of injury and falls with supportive stable positioning




Comfort Tension Seating[®] increases support as well as sitting tolerance

Enhanced positioning with refined features

UPHOLSTERY OPTIONS



FRAME OPTIONS



OPTIONAL ACCESSORIES



Foot Box



Lateral Support



Tray



O² Holder



Half Tray



Padded Thigh Strap

SPECIFICATIONS

Degrees of Tilt	Up to 38°
Recline	Up to 90° (Trendelenburg capable)
Weight Capacity	up to 350 lbs. (114 kg)
Seat Widths	16", 18", 20", 22", 24", 26", 28" (40cm, 45cm, 50cm, 55cm, 61cm, 66cm, 71cm)
Seat Depths	15.5", 17", 18.5", 20" (39cm, 43cm, 45cm, 58cm)
Floor-to-Seat Heights	Adjustable 17", 18", 19", 20", 21", 22" (43cm - 55cm)
Wheels	5" Front Casters: Total Lock (12.5cm) 5" Rear Casters: Directional or Total Lock (20cm)
Lifetime Warranty Frame	10 years on 16 Gauge tubular steel frame
Pressure-Reducing Seat & Backrest	Integrated with Comfort Tension Seating [®]
Standard Padding	Neckrest, shoulder bolsters, upper and lower side pads, seat and back pad, calf and sole pad



**ProCare
Medical**

ProCare Medical

759 Flory Mill Road, Lancaster PA 17601

www.procare-medical.com

Ph: 855-528-0421 | Fax: 877-528-0421