

# Encore

*Comfort, Functionality, and Style*

**BRODA®**



*New frame design, including adjustable wings.*



*Enhanced mobility and tilt-in-space positioning gives users freedom.*



*Pedal Rocker calms and creates contentment.*



*Comfort Tension Seating® increases support and sitting tolerance.*

## ***Patient, Facility, and Caregiver Benefits***

- *Increases quality of life and engagement with surroundings.*
- *Allows for easy self-mobility and safer patient maneuverability.*
- *Reduces risk of injury and falls with supportive stable positioning.*

### The Encore of Comfort, Functionality, and Style.

#### SPECIFICATIONS


Tilt - Infinitely Adjustable	<b>Rocker:</b> up to 30° (combined with rocking) <b>Chair:</b> up to 25°
Weight Capacity	<b>Rocker:</b> up to 250 lbs. / 158 kg <b>Chair:</b> up to 350 lbs. / 114 kg
Seat Widths Available	<b>Rocker:</b> 16", 18", 20", 22" (40cm, 45cm, 50cm, 55cm) <b>Chair:</b> 16", 18", 20", 22", 24", 26" 28" (40cm, 45cm, 50cm, 55cm, 61cm, 66cm, 71cm)
Seat Depths Available	<b>Rocker:</b> 15.5", 17", 18.5", 20" (39cm, 43cm, 46cm, 50cm) <b>Chair:</b> 17", 18.5", 20" (43cm, 46cm, 50cm)
Back Lengths Available	22", 32", 33.5", 35" (55 cm, 81cm, 85cm, 89cm)
Floor-to-Seat Heights Available	<b>Rocker:</b> Adjustable 14" - 19" <b>Chair:</b> Adjustable 12.5" - 17.5"
Wheels - Casters	5" front total lock (12.5cm) 5" fixed (12.5cm)
Wheels - Mag Wheels	20" or 22" front mount (50 cm or 55 cm)
Lifetime Warranty Frame	10 years on 16 Gauge tubular steel frame
Pressure-Reducing Seat & Backrest	Integrated with Comfort Tension Seating®, Solid Seat Pan Available
Standard Padding*	Headrest cushion, shoulder bolsters, lower side pads, seat and back pad with cover

\*22" back option excludes headrest cushion and shoulder bolsters.

#### UPHOLSTERY OPTIONS



#### FRAME OPTIONS



Shown with full padding and adjustable wings

Additional accessories, such as foot supports, specialty padding packages, trays, and positioning belts are available.

