

Where Comfort is Essential



Pedal Chair - Fully Padded



A Superior Mobility Chair!

- Broda Comfort Tension Seating™ helps **reduce pressure and heat and moisture build up**, creating long term sitting comfort.
- 18° of infinitely adjustable tilt in space cradles the body, **preventing slumping and sliding without the use of restraints and is exceptionally comfortable.**
- Front pivot seat tilt allows for a **proper foot on floor position** for foot propelling while tilted- **Encourages safe mobility and independence.**
- Adjustable seat to floor height from 12.5" to 17.5" **accommodates both short and tall residents.**
- Easily adjustable push pin arm height adjustment enables **proper arm positioning** for various residents.
- 16 gauge tubular steel frame ensures **durability** and helps **reduce chair maintenance and replacement cycles.**

Broda Comfort Tension Seating™



- Reduces pressure, thus helping to reduce skin breakdown
- Creates air flow for heat and moisture exchange
- Memory retention allows strapping to return to its original shape within seconds

Where Comfort is Essential

Excellent for Fall Prevention!



The ergonomic design maintains a proper "Foot-On-Floor" position when using the seat tilt, helping to prevent falls

Standard Features

- Broda Comfort Tension Seating[™]
- 16 gauge tubular steel powder coated frame
- Up to 18° of infinitely adjustable seat tilt
- 18" or 20" standard seat widths
- 17" seat depth
- 12.5" to 17.5" adjustable seat height
- 32" back length
- Neckrest cushion
- Lower side pads
- Adjustable arm height
- 5" Total lock front casters, 5" fixed rear casters
- 350 lb. weight capacity

Options

- Full padding upgrade
- Mag wheels (rear or front mounted)
- Swing away and/or elevating swing away footrests
- ABS solid tray, clear tray, padded tray cover
- Vinyl covered seat and back pad
- IV Pole
- O2 Holder
- Lateral supports, positioning cushions
- Custom and Bariatric sizes available



Vinyl covered seat and back pad, and clear tray – Clear tray allows agitated residents to see their feet



Optional 20" mag wheels
Excellent for residents who can propel with their arms



**ProCare
Medical**

The Pedal Chair provides superior comfort, and versatility making it ideal for various settings:

- Long Term Care
- Acute Care
- Hospitals
- Home Care