

Revive

Shower Commode

HYGIENE
BRODA[®]



Anterior and posterior tilt aid in safe patient handling



Tilting and reclining permits safe and thorough bathing



Enhances patient comfort and support with adjustable lower leg support



Removable arm supports allow for safe and easy patient transfers

SPECIFICATIONS

Degrees of Tilt	Up to 12° posterior Up to 10° anterior
Degrees of Recline	Up to 67°
Seat Widths	18", 22" (45mm, 55mm)
Seat Depths	18" (45mm)
Weight Capacity	400 lbs. (181kg)
Lifetime Warranty Frame	10 years on 16 gauge tubular stainless steel white frame



Concave back option accommodates and supports gluteal shelf



Dynamic tilting and reclining permits safe and thorough bathing



Anterior and posterior tilt allows for easy user positioning



Rust-resistant frame and components are durable and long-lasting



Shown with optional head support

SPECIFICATIONS

Degrees of Tilt	Up to 14° posterior Up to 8° anterior
Degrees of Recline	Up to 65°
Seat Widths	27", 33" (68mm, 83mm)
Seat Depths	21", 23" (53mm, 58mm)
Weight Capacity	800 lbs. (362kg)
Lifetime Warranty Frame	10 years on 16 gauge tubular stainless steel white frame



**ProCare
Medical**

ProCare Medical, 759 Flory Mill Road, Lancaster PA 17601
Ph: 855-528-0421 | www.procare-medical.com | Fax: 877-528-0421