



ROPOX VARIO CHANGING BED

SPACESAVING CHANGING BED

The Ropox VARIO changing bed is a full size height adjustable unit mounted to a wall.

The height is adjustable from 12" -39.5" using the hand control. The control box plugs into a standard electrical outlet.

When the changing bed is not in use, the tabletop can be folded up to a vertical position so that it takes up minimal space in the room.

Folding up and down the table top takes place with easy movement and pneumatic cylinder assist.

Ropox Vario changing bed is available in five models including and without a folding side bed guard and mattress adapted to the changing table.

The mattress has a cleaning-friendly surface that makes it easy to clean after use.

Also available is a shower and changing table, see separate literature.



ROPOX™






BETTER WAYS TO BETTER DAYS






- A PART OF AddLife®

ROPOX VARIO CHANGING BED

Height adjustment: 27-1/2" (11-3/4"-39-1/4")
 User weight: 440 lb DS/EN 12182:2012

ROPOX VARIO CHANGING BED

VARIO WITHOUT BED GUARD RAIL				
				
Sized to fit 5' (60") space				
47-1/4" x 27-1/2" 40-30604	55" x 27-1/2" 40-30606	63" x 27-1/2" 40-30608	70-3/4" x 27-1/2" 40-30610	74-3/4" x 27-1/2" 40-30611


VARIO WITH BED GUARD RAIL				
				
Sized to fit 5' (60") space				
47-1/4" x 27-1/2" 40-30804	55" x 27-1/2" 40-30806	63" x 27-1/2" 40-30808	70-3/4" x 27-1/2" 40-30810	74-3/4" x 27-1/2" 40-30811

* (3) sizes stocked, the others are available with lead time

ACCESSORY

Cushion to the bed:

- 40-25782 for 47-1/4" / 55"
- 40-25783 for 63" / 70-3/4" / 74-3/4"



Mattress tabletop (thickness 1-1/4")

- 40-25760 for 47-1/4"
- 40-25795 for 55"
- 40-25796 for 63"
- 40-25764 for 70-3/4"
- 40-25765 for 74-3/4"



ProCare is the exclusive distributor of Ropox products in the United States. Contact our office or your local dealer today for more information.

www.procare-innovations.com 855-528-0421

