



INTRODUCING THE MOTORIZED

IMMO WEIGHT CART SYSTEM

A compact solution utilized by hospitals systems, biomedical departments and leading ceiling lift companies around the world to load test lift systems.

DESIGNED FOR HEALTHCARE

- 3 models to choose from
- Designed for testing of new installations or regular periodic maintenance
- UL and CE approved

STANDARD FEATURES

- Power drive to eliminate any manual pushing
- Selectable weight to match capacity of lift system
- Battery and charger included



EFFICIENCY IN OPERATIONS

Discover the innovative method to weight test lift systems in any room or setting. The cart system can be utilized by biomedical and facility technicians immediately when a room is vacant or available. Also suitable for mobile lifters and can be transported easily in a van.



SAFETY MEASURES

With increased utilization of lift systems, the testing of these systems can cause musculoskeletal injuries to technical staff if manual weight plates and carts are utilized. Provide a safe solution and comply with NIOSH AND OSHA recommendations for workplace ergonomics.



Integrates with any brand of lift equipment

“The IMMO cart system has been a great improvement in the safety and process of inspecting our patient lifts. Highly recommend it for routine maintenance and eliminating risk of lifting 45 lb plates and pushing a cart.”





Click or scan to watch video

	Model 255	Model 350	Model 550
Capacity	561 lb 255 kg	770 lb 350 kg	1,210 lb 550 kg
Dimensions	Width: 21.5" Depth: 26.4" Height: 37.8"	Width: 31.5" Depth: 26.4" Height: 47.2"	Width: 31.5" Depth: 26.4" Height: 47.2"
Motor and Drive	450W	1200W	1200W
Anti-Tip	Yes	Yes	Yes
Application and Notes	Smallest size unit, ideal for mobile lifts and basic ceiling lifts	These two models have same footprint, the 550 model has more weight plates and upgraded chassis, motor, and drive system	



ACCESSORIES:

- #101557 Extra Battery Pack
- #102255 Wall Charger Station and Battery Pack
- #100613 Strap and hardware kit
- #101501 Adjustment Weight kit (5 kgs and 1 kgs)

ABOUT US:

ProCare is the exclusive provider of the IMMO cart systems for USA and Canada.

IMMO is a global producer of powered carts and trolley systems since 1967.

CONTACT US

Phone: 855-528-0421