




ProCare Innovations

transforming the shower experience
protecting the caregiver.



Includes DURA[®] care antimicrobial coating*
Best in class 5 year warranty



Healthcare facilities are increasingly utilizing showering for users' personal hygiene care. Showering a user in a standard chair will cause the caregiver to work in unnatural postures and suffer static back strain. The best way to overcome this situation is to introduce the proper equipment that significantly **increases the percentage of time that the caregiver spends working in a healthy position**. Lopital offers comfort and safety with an extensive line of shower chairs and trolleys.

The various models of the adjustable and hi-low shower chairs provide the user maximum safety and comfort. A **single caregiver can** perform dressing and undressing, showering, toileting, foot care, and hair washing in a proper ergonomic working condition. Lopital shower chairs are designed to be utilized across the spectrum of mobility function, from active to semi and fully dependent.

The **Lotus** is an innovative active stand aid to assist users from seated to standing positions and transfers. The **Marina** shower trolley provides an ergonomic solution for showering dependent users with no ability to support themselves. Users are showered in a supine position on a comfortable soft mattress. Transfers to the shower trolley from the bed can be performed without manual lifting using a ceiling lift, mobile lift, or horizontal transfer device.

Product Details

Product no.	Product description	Adjustable hydraulic/ electric or fixed	Mobile	Height adjustable	Headrest	Tilt mechanism	Shower	Toilet	Closed cell foam	Tech. lifespan min. 10 years*	Armrests with front locking	Can be wheeled over toilet	Bariatric
5100 5199	Flexo Shower-toilet chair	hydraulic	•	•	optional	*	•	•	•	•	•	•	
5100 5500	Elexo Shower-toilet chair	electric	•	•	optional	*	•	•	•	•	•	•	
5100 5300	Elexo XXL Shower-toilet chair	electric	•	•		*	•	•	•	•		•	•
5100 5600	Reflex Shower-toilet chair	electric	•	•	•	•	•	•	•	•	•	•	
5100 5100	Nemo Tilt in Space S/T	hydraulic	•		•	•	•	•	•	•	optional	•	
5100 5705	Tango XL Shower-toilet chair	fixed	•				•	•	•	•		•	•
5100 5700	Tango Shower-toilet chair	fixed	•				•	•	•	•	optional	•	
5100 5605	Reflex Shower-toilet chair with foldable footrests	electric	•	•	•	•	•	•	•	•	•	•	
5100 5800	Tango XXL	fixed	•				•	•	•	•		•	•
6100 2000/2010	Marina Basic or Deluxe	hydraulic	•	•	•	•	•			•			•
6100 2010/2060	Marina Basic or Deluxe	electric	•	•	•	•	•			•			•

*Tilts automatically with height adjustment.



5100 5100
Nemo Tilt in Space Shower-Toilet Chair
Capacity: 330 Lbs.
Seat Width: 18"



5100 5199
Flexo Shower-toilet chair
Capacity: 330 Lbs.
Seat Width: 18"



5100 5500
Elexo Shower-toilet chair
Capacity: 400 Lbs.
Seat Width: 18"



5100 5300
Elexo XXL Shower-toilet chair
Capacity: 485 Lbs.
Seat Width: 28"



5100 5600
Reflex Shower-toilet chair
Capacity: 300 Lbs.
Seat Width: 18"



5100 5605
Reflex Shower-toilet chair with foldable footrests
Capacity: 300 Lbs.
Seat Width: 18"



5100 5700
Tango Shower-toilet chair
Capacity: 400 Lbs.
Seat Width: 18"



5100 5705
Tango XL Shower-toilet chair
Capacity: 550 Lbs.
Seat Width: 22"



5100 5800
Tango XXL Shower-toilet chair
Capacity: 800 Lbs.
Seat Width: 28"



6100 2000/2010 Hydraulic
6100 2050/2060 Electric
Marina Shower trolley
Capacity: 500 lbs

* Technical life span is only guaranteed with correct use and maintenance in accordance with the supplier's instructions.



Reflex XL

5100 5950

- 500lb weight capacity
- 30 degree tilt-in-space
- 22" seat width
- Adjustable backrest, headrest, footrest, and armrest

Reflex XXL

5100 5900

- 600lb weight capacity
- 30 degree tilt-in-space
- 28" seat width
- Adjustable backrest, headrest, footrest, and armrest



Tango Series



Tango XL

5100 5705N

- 550lb weight capacity
- 22" seat width

Tango XXL

5100 5800



- 800lb weight capacity
- 28" seat width

Lotus Active Stand Aid

Product no.	Product description	Material	Lifting capacity lbs	Length x width x height	Width legs	Mechanical base spreading	Weight lbs
5900 8200	LOTUS	Steel	440 lbs	39.4" x 26.8" x 43.3"	26.8"	26.8" - 37.8"	77
5900 8250	LOTUS COMPACT	Steel	440 lbs	39.4" x 26.8" x 39.4"	26.8"	26.8" - 37.8"	77
5900 8270	LOTUS XL	Steel	500 lbs	39.4" x 28" x 43.3"	26.8"	26.8" - 37.8"	77



Cross-infection issues can be reduced or eliminated by equipment properly designed to ease the caregiver's task ensuring a thorough shower is given. Equally as important to reducing cross-infection is equipment designed to be easily disinfected post shower or hygiene activity. Lopital products are designed to meet that requirement.

The user enjoys comfort and dignity and the caregiver benefits from the improved safety and ergonomics of the Lopital product. This leads to a higher quality of care and a superior showering and hygiene experience for the user.

Reduce Cross Contamination

Quality of Life

Positive Clinical Outcomes

Safety and Ergonomics



Everybody involved in the caregiving venue strives for a positive outcome for their patient or resident. When you bring all the factors identified above together with the proper tools, Lopital empowers you to improve outcomes and enhance your overall quality of care.

Your showering and hygiene equipment should reduce lost time, reduce lost time injuries, and the cost of care for your facility. Lopital integrates into your safe patient handling program through the thoughtful ergonomic design and function of our products.

A complete range of solutions for improving outcomes



**ProCare
Innovations**

www.procare-innovations.com

759 Flory Mill Road Lancaster,
PA 17601 Phone: 855-528-0421
Fax: 877-528-0421

FSS Contract #36F79720D0131